

SENTINELTM AccuTOUCH

Powder-free Latex Surgical Gloves

Unique **PRO-FitTM design** ensures an accurate fit, maximum dexterity

Superior grip performance due to **anti-slip smooth finish**

Biomechanically superior design with unconstrained thumb, results in lesser hand fatigue

Enhanced tactile sensitivity due to thinner design

Hydrogel coating enhances wet donning ability

AQL* 0.65 ensures reliable and higher levels of protection

SENTINEL™ ACCUTOUCH

Powder-free Latex Surgical Gloves

***AQL = Acceptance Quality Limit.**

In this case, it refers to the air leak test for detection of holes, which is a critical standard for assuring surgeons that they are being protected from holes, rips or very weak areas. Our AQL 0.65 reflects a tighter quality assurance level, giving the surgeon an assured higher level of protection.

Other gloves are at least 3.5 times as likely to fail than Sentinel™ brand gloves.

Elite Surgical Sentinel™ AccuTouch surgical gloves feature a unique PRO-Fit™ design - a superior 3-layer ultra-thin design, featuring two layers of latex with a hydrogel coating - that ensures an accurate and natural fit, provides maximum dexterity and results in lesser hand fatigue. The ultra thin design also allows for enhanced tactile sensitivity.

The anti-slip smooth finish of Sentinel™ AccuTouch offers a remarkably superior grip performance enabling natural hand movement. The low protein content ensures low levels of allergic contact dermatitis.

Free the Thumb

The unconstrained thumb design with straight fingers places the thumb away from, or in front of the finger plane, enabling ease of hand movements during surgery.

Product Description Latex Wet Donning Sterile Surgical Gloves, powder-free, smooth surface, straight fingers, beaded cuff, size embossed.

Product Specifications

REF	ES-SAT							
Size (mm)	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
Palm Width (mm)	70±6	76±6	83±6	89±6	95±6	102±6	108±6	114±6
Overall length								
mm	310 ± 10							
inch	12.2 ± 0.4							
Single Thickness	Finger		Palm		Cuff			
mm	0.19 ± 0.02		0.18 ± 0.02		0.15 ± 0.01			
mil	7.48 ± 1.19		7.09 ± 0.79		5.91 ± 0.39			
Sterilisation	Gamma irradiation							
Quality Standards	Conforms to ASTM D3577, ASTM D5712, EN455, ISO 10282 and AS/NZS 4179 Certified with CE, FDA							
AQL	Pinhole AQL 0.65 as per EN455-1 / ASTM D5151							
Packaging	50 pairs/box; 6 boxes/carton							
Shelf life	Three (3) years from date of manufacture							



Elite Surgical Pty Ltd
41 Parriwi Road, Mosman, NSW 2088
Australia

E-mail: info@elitesurgical.com.au
Website: www.elitesurgical.com.au