


SUCTION & IRRIGATION

Quality healthcare products for hassle-free procedures



Features

Tubing 215 cm, Latex free
PVC (DEHP free)

PRODUCT CODE	ES080601 2-way, 1 Chamber
Non-vented Spike	2
Drip Chamber	1
Click Clamp	3
Y-Connector	1
Silicone Tube Adaptor	1

Dual perforators with click clamps, trauma pump, click clamp, stepped connector with silicone adaptor.

The silicone adaptor is attached to the arthroscope. The trauma pump chamber and click clamp allow for easy control of flow rate.




Dual perforators with click clamps, two trauma pumps, click clamps, stepped connector with silicone adaptor.

The silicone adaptor is attached to the arthroscope. The dual chambers facilitate better, faster and more regulated irrigation; click clamps and allow for easy control of flow rate.

Features


Tubing 215 cm, Latex free
PVC (DEHP free)

PRODUCT CODE	ES080602 2-way, 2 Chamber
Non-vented Spike	2
Drip Chamber	2
Click Clamp	3
Y-Connector	1
Silicone Tube Adaptor	1



SUCTION & IRRIGATION

Quality healthcare products for hassle-free procedures



ARTHROSCOPY IRRIGATION SET

PRODUCT CODE	ES080603
Non-vented Spike	Four spikes
Drip Chamber	Single Pump
Click Clamp	260 cm
Y-Connector	Double Packaging
Silicone Tube Adaptor	EO sterile