

# SENTINEL™ ORTHO

Powder-free Latex Surgical Gloves

Unique **PRO-Fit™ design** ensures an accurate fit, maximum dexterity and lesser hand fatigue

**Superior grip** performance due to **textured fingertips**

**Extra thickness** for ultimate protection during orthopaedic procedure

**Hydrogel coating** enhances wet donning ability

**Brown color** reduces glare from operating lights

# SENTINEL™ ORTHO

Powder-free Latex Surgical Gloves

**\*AQL = Acceptance Quality Limit.**

In this case, it refers to the air leak test for detection of holes, which is a critical standard for assuring surgeons that they are being protected from holes, rips or very weak areas. Our AQL 0.65 reflects a tighter quality assurance level, giving the surgeon an assured higher level of protection.

Other gloves are at least 3.5 times as likely to fail than Sentinel™ brand gloves.

Elite Surgical Sentinel™ ORTHO surgical gloves feature a unique PRO-Fit™ design - a superior 3-layer design, featuring two layers of latex with a hydrogel coating - that ensures an accurate and natural fit, provides maximum dexterity and results in lesser hand fatigue. Extra thickness provides extra extra protection during rough use in orthopaedic procedure

The textured fingertips offers a remarkably superior grip performance enabling natural hand movement. The very low protein content - less than 50µg/dm<sup>2</sup> - ensures low levels of allergic contact dermatitis (Type IV allergy).

Made from select natural latex and offered in Brown colour, Sentinel™ ORTHO gloves offer excellent wet donning ability, and provide optimum protection and a comfort fit throughout the surgical procedure.

**Product Description** Latex Wet Donning Sterile Surgical Gloves, powder-free, textured fingertips Curved fingers, non-beaded cuff, size embossed.

## Product Specifications

REF	ES-SOP							
Size (mm)	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
Palm Width (mm)	70±6	76±6	83±6	89±6	95±6	102±6	108±6	114±6
<b>Overall length</b>								
mm	310 ± 5							
inch	12.2 ± 0.2							
<b>Single Thickness</b>	<b>Finger</b>		<b>Palm</b>			<b>Cuff</b>		
mm	0.34 ± 0.03		0.33 ± 0.02			0.30 ± 0.01		
mil	13.39 ± 1.18		12.99 ± 1.18			11.89 ± 1.18		
<b>Sterilisation</b>	25 kGy Gamma irradiation							
<b>Quality Standards</b>	Conforms to ASTM D 3577, ASTM D 5712, EN455, ISO 10282 and AS/NZS 4179 Certified with CE, FDA , ASTM F 1671							
<b>AQL</b>	Pinhole AQL 1.5, 1.0, 0.65 (optional)							
<b>Packaging</b>	50 pairs/box; 6 boxes/carton							
<b>Shelf life</b>	Three (3) years from date of manufacture							



Elite Surgical Pty Ltd  
Unit 49B, 2 Slough Avenue,  
Silverwater NSW2128, Australia

**E-mail:** info@elitesurgical.com.au  
**Website:** www.elitesurgical.com.au